

# Preclinical irradiations at conventional and very high dose rates (FLASH-RT) on small animals at the ARRONAX research cyclotron: repositioning, margins and field sizes

Manon Evin, Charbel Koumeir, Gaëlle Saade, Mathieu Chocry, Ferid Haddad, Vincent Potiron, Noël Servagent, Daphnée Villoing, Stéphane Supiot, Vincent Métivier, et al.

#### ▶ To cite this version:

Manon Evin, Charbel Koumeir, Gaëlle Saade, Mathieu Chocry, Ferid Haddad, et al.. Preclinical irradiations at conventional and very high dose rates (FLASH-RT) on small animals at the ARRONAX research cyclotron: repositioning, margins and field sizes. 61èmes Journées Scientifiques de la Société Française de Physique Médicale (SFPM), Jun 2023, Nancy, France. . hal-04145654

# HAL Id: hal-04145654 https://imt-atlantique.hal.science/hal-04145654v1

Submitted on 29 Jun 2023

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers.

L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.









# ronax lantes

#### Preclinical irradiations at conventional and very Sfpm high dose rates (FLASH-RT) on small animals at the ARRONAX research cyclotron: repositioning, margins and field sizes



Manon Evin a\*, Charbel Koumeir ac, Gaëlle Saade d, Mathieu Chocry bd, Ferid Haddad ac, Vincent Potiron bd, Noël Servagent a, Stéphane Supiot bd, Daphnée Villoing b, Vincent Métivier a, Gregory Delpon ab and Sophie Chiavassa ab

a Nantes Université, IMT Atlantique, CNRS/IN2P3, SUBATECH, F-44000 Nantes, France, b Institut de Cancérologie de l'Ouest, Saint-Herblain, France, GIP ARRONAX, Saint-Herblain, France, d Nantes Université, CNRS, US2B, UMR 6286, F-44000 Nantes, France

#### Context

In the framework of translational research, it is envisioned to irradiate mice on radiobiology platforms, as the one of the ARRONAX cyclotron. In general, these facilities do not feature an integrated imaging system, which makes irradiation of organs without external anatomical landmarks challenging. Moreover, the field size cannot be tailored to each animal but must be set for a cohort.

This study aims to address the parameters required to determine the optimal field size for a selected cohort, which can be expressed as:

 ${\it Optimal Field size (cohort) = Field size (cohort) + \; Margins \; }_{\it Targeting} + \; Margins \; {\it Respiratory}$ 

### (1) Targets dimensions and inter-animal variability

The volumes of the boxes containing the brain, the lungs and the intestines were measured on 3x3 distinct tomographic images (XRAD225Cx,  $\mu CBCT$ , 40kVp, 0.2mm isotropic voxel size) for 11 mice (C57BL/6 - 5 females, 6 males) and compared to the animal's weight, as shown in Fig. 1. The relative uncertainty on the volume is the quadratic sum of the uncertainty of each length, corresponding to the size of two pixels.

- Brain volume was found to be independent of the weight
- Unique field size
- Lungs and intestinal volumes are correlated with the mouse weight
- Selection of a cohort on the weight criterion enables a unique field size




Figure 1. For the three locations, volume containing the organ with respect to the animal weight.

## Immobilization and targeting system

The system, shown in Fig. 2, consists of a patient-like 3D-printed bed and harmless contentions, together with graduations as external markers for the lasers.

The targeting accuracy was assessed by resetting 10 times an anaesthetized C57BL/6 mouse for the brain, lungs and intestines. Each time, a  $\mu CBCT$  (40kVp, 0.2mm isotropic voxel size) was acquired using the XRAD225Cx. The TPS µRaystation® was used for organ contouring on the images and to perform rigid registrations, without rotation, with respect to the reference image.

The maximum offsets in the three directions were recorded in the Table 1. To ensure the target coverage, equivalent margins should be applied for each organ. An example of the histogram of deviations obtained for a site is also

	Right/Left	Longitudinal	Antero-Posterior
Brain	1,4	0,6	0,9
Lungs	2,0	3, 4	0, 4
Intestines	2, 3	3,0	0, 7

Table 1 – Maximum offsets with respect to the reference image (mm)



Figure 2. Photograph of the mouse holding and targeting system




Figure 3. Histogram of brain offsets measured on the 10 images.

## Respiratory motion -

In addition to the accuracy of the targeting, margins should also take into account potential organ motions such as respiratory motion. The fluoroscopic mode of the XRAD225Cx (40kVp, 0.1 mm voxel size) was used to acquire 10 respiratory cycles for a mouse. The wider cycle is displayed in Fig. 4.

Measured maximum amplitudes (full insipiration) relatives to full expiration position (Fig. 5.) are:

- · 1.7 mm in the longitudinal direction (diaphragm)
- 0.6 mm in the antero-posterior direction.




Figure 4. Maximal amplitude cycle found with respect to time. Interpolation with Akima method.



Figure 5. Fluoroscopic mouse views

Moreover, the respiration duration, which lasts 2.3 s in total, is long in comparison with FLASH-RT irradiation time (few hundred milliseconds) [1]. Then, the above amplitudes must be added to the margins to ensure the lungs coverage, where a loss could occur for 30 % of the irradiation time in conventional dose rate and up to 100 % of the time in FLASH-RT.

#### Conclusion and perspectives

- The margins obtained and the study of inter-animal variability guarantee optimal field sizes for irradiation of the brain, lungs and intestine.
- The methodology developed ensures accurate and reproducible positioning for targeted preclinical irradiations at ARRONAX.

Following steps will involve Monte Carlo GATE simulations to:

- Evaluate of the dosimetric impact of the respiratory cycle using the 4D digital mouse fantom (MOBY) [2],
- · Verify the organ coverage, as well as distal and lateral dose homogeneity, on the cohort with adapted field size and margins found in this study.